

Why to use ELVeS® Radial® Swift Fibers?

- ___ Just 5Fr Introducer Size
- __ or Easy access via 14GA catheter

Well known FUSION® safety concept

The swift fibers are designed to complete our product range. You decide upon access by direct puncture via 14GA canula or using an ordinary 5Fr introducer set (instead of 6Fr).

Adapt flexibly to the particular vein in terms of energy, power and approach.



Placing the 14GA catheter



Entering the vein



Starting laser treatment

Fibers*

REF	Product	PU **	length [m]	ø fiber tip [mm]	CI***
503100140	ELVeS® Radial® swift, IC	10	2.5	1.5	5Fr/14GA
503100142	ELVeS® Radial® swift 2ring, IC	10	2.5	1.5	5Fr/14GA

Kits*

503100141	ELVeS® Radial® swift Kit, IC	10	2.5	1.5	5Fr/14GA
503100143	ELVeS® Radial® swift 2ring Kit, IC	10	2.5	1.5	5Fr/14GA

^{*} Only to be used with 1470 nm ** Packaging unit *** Compatible introducer

Imprint
biolitec AG
Untere Viaduktgasse 6/9
A-1030 Vienna
Phone: +43 1 3619 909 50
info@biolitec.de
www.biolitec.com

All fibers are free of latex and DEHP.

Our fibers are single use products
(unless otherwise indicated) delivered
sterile for immediate use.
biolitec®, FUSION® and ELVeS® Radial®
are registered trademarks owned
by biolitec.



