

LaEvita

Regain your femininity

The laser treatment of genito-urinary syndrome of the menopause

biolitec® has the complete solution to the main disorders affecting the female reproductive system.



Why LaEvita?

LaEvita is the new virtually pain-free solution to the treatment of **vaginal atrophy** and **mild stress urinary incontinence**. Unlike other CO₂-lasers, this is a **NON-ablative and NON-destructive tissue technique** that is virtually painless and well tolerated by patients.

How does it work?

Thanks to the special applicator, the LEONARDO® Dual laser emits a gentle, painless light that stimulates collagen production and restores vaginal elasticity and hydration, and improves vascularization of the vaginal mucosa. Each session lasts about 10 minutes and is performed without anesthesia. The treatment includes 3 sessions at 4-week intervals. In more severe cases, 1 additional session may be necessary. The procedure is performed on an outpatient basis, is virtually painless for the woman, and does not require hospitalization.

What is required?

The technique requires the use of the LEONARDO® Dual laser, the use of special reusable LaEvita handpieces and dedicated kits for the treatment of vaginal atrophy or urinary incontinence

Advantages

- It significantly reduces vaginal dryness and pain without side effects
- It helps maintain a fulfilling intimate life
- Fully outpatient procedure
- No anesthesia required
- No side effects
- Effective
- Almost no pain
- No downtime



LaEvita Handpieces*

REF	Product	Outer ø [mm]	Inner ø [mm]	Length ø [mm]
400500400	LaEvita Small	24	2,20	100
400500410	LaEvita Medium	24	2,20	135
400500420	LaEvita Large	28	2,20	135
400500430	LaEvita Handpieces Set - (LaEvita small, medium, large)			

*The Vaginal Handpiece is a reusable for use in gynecological vaginal treatments.



LEONARDO® DUAL 45
The complete solution
for gynecology and more

Kits

REF	Product	Components	Quantity
503200610	Vaginal Atrophy LaEvita Kit, IC	Coupling Fiber, IC	1
		Side Fiber Tip	10
503200620	Urinary Incontinence LaEvita Kit, IC	Coupling Fiber, IC	1
		Radial Fiber Tip	10
		Side Fiber Tip	10



biolitec biomedical technology GmbH

Phone: +49 3641 519 53 0

Otto-Schott-Str. 15

07745 Jena, Germany

Contact us

